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
● Scheduled Use

The heating plate IKATHERM HP 30 digital is a universally applicable laboratory heating instrument. Its temperature can be adjusted infinitely variable from ambient temperature up to +400°C and it is controlled by a cast integral thermocouple. The instrument has a provision for connecting an external sensor.

● Unpacking

Please unpack the instrument carefully and look out for any damage. It is important to recognize a transportation damage already when unpacking. If this is the case, an immediate factual report is required (mail, rail or forwarding agency).

● Setting-up the Instrument

Attention has to be paid to the ambient conditions listed in technical data (temperature and humidity).

Safety Warnings:

The instrument must not be operated in rooms with explosion hazards.

The instrument has to be placed upon a fire-proof and/or not burnable surface.

The mains cable may not get in contact with the heating plate!


The heating plate heats up to 400°C. For this reason it is dangerous to touch the heating plate or the guard plate.

● Starting-up

Before starting up the instrument, please check whether the voltage indicated on the ID-plate of the instrument corresponds with the mains voltage at hand. If this is not the case, the instrument must not be operated. The instrument has to be plugged into an earthed socket.

The heating plate is switched on at the rear with the main switch.

● Operating and Function Elements

After the instrument has been switched on via the main switch, always status "%", "process value" and "intern sensor" are displayed via the red LEDs above the respective keys.

1) Display	Indication of "real" or "rated"-temp.
2) Process value	Indication of "real" temperature on the display, measured by the built-in or by the external temperature sensor.
3) Set value	Display of the "rated" temperature on the display, adjusted by means of the control dial (11).
4) 0 %	Corresponds with "stand-by"-operation, but the heating plate is not heated-up.
5) 25 %	This key switches on 25 % of the heating power.
6) 50 %	This key switches on 50 % of the heating power.
7) 100 %	This key switches on 100 % of the heating power.
8) Set Temp	The requested "rated" temperature is set with the control dial (11).
9) Intern Sensor	In the mode "process value", the "real" temperature of the heating plate is displayed. The control is carried out via the built-in temperature sensor.
10) Extern Sensor	This function is only possible with the external sensor HP 30.1 that is connected on the rear of the instrument to a 6-polar bushing. In the mode "process value", the temperature measured on the sensor is displayed. The temperature control is carried out via the external measuring sensor. By means of the built-in internal temperature sensor the heating plate temperature is limited to the maximum adjustable value +10°C.
11) Control dial	The requested "rated" temperature is set via the control dial.

● Heating

The temperature of the heating plate is kept constant by an electronic controller. The temperature sensor required for this, is built-into the heating plate.

The instrument has a cast iron heating plate with a heating power of 2500 watts (at 100 %).

The requested heating power can be selected.

In order to prevent an overshooting of the temperature, it is recommended to reduce the heating power to 25 %, if for example small vessels are used.

The heating process can be interrupted any time by the "0 %"-key.

The temperature of the heating plate is limited to 300°C when using an external temperature sensor.

In case of a defect of the installed or of the external temperature measuring sensor, the display shows "452".

● Maintenance and Cleaning

The instrument requires no maintenance. In case of repair, it may only be opened by a qualified technician.

The outside of the instrument may be cleaned with:

- a) water with a tenside containing detergent
- b) isopropyl alcohol (e.g. if very soiled)

Please submit only instruments for repair that have been cleaned and are free from substances dangerous to health.

When ordering spare parts, please quote type of instrument as well as serial number.

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● Guarantee

You have acquired an Original IKA - Laboratory Instrument that meets highest requirements in technology and quality.

In accordance with the IKA - terms of warranty, this instrument is guaranteed for a period of 12 months.

In the event of a claim under guarantee, please contact your laboratory dealer. You can also send the instrument directly to our Works. Please attach freight invoice and state your complaints. Freight charges are at your expense.

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