IKATRON® ETS-D6

<i>BETRIEBSANLEITUNG</i>	<i>DE</i>	<i>7</i>
<i>OPERATING INSTRUCTIONS</i>	<i>EN</i>	<i>16</i>
<i>MODE D'EMPLOI</i>	<i>FR</i>	<i>25</i>
<i>ADVERTENCIAS DE SEGURIDAD</i>	<i>ES</i>	<i>34</i>
<i>VEILIGHEIDSIINSTRUCTIES</i>	<i>NL</i>	<i>35</i>
<i>NORME DI SICUREZZA</i>	<i>IT</i>	<i>37</i>
<i>SÄKERHETSANVISNINGAR</i>	<i>SV</i>	<i>38</i>
<i>SIKKERHEDSHENVISNINGER</i>	<i>DA</i>	<i>40</i>
<i>SIKKERHETSANVISNINGER</i>	<i>NO</i>	<i>41</i>
<i>TURVALLISUUSOHJEET</i>	<i>FI</i>	<i>43</i>
<i>INSTRUÇÕES DE SERVIÇO</i>	<i>PT</i>	<i>44</i>
<i>WSKAZÓWKI BEZPIECZEŃSTWA</i>	<i>PL</i>	<i>46</i>
<i>BEZPEČNOSTNÍ POKYNY</i>	<i>CS</i>	<i>47</i>
<i>BIZTONSÁGI UTASÍTÁSOK</i>	<i>HU</i>	<i>49</i>
<i>VARNOSTNA OPOZORILA</i>	<i>SL</i>	<i>50</i>
<i>BEZPEČNOSTNÉ POKYNY</i>	<i>SK</i>	<i>52</i>
<i>OHUTUSJUHISED</i>	<i>ET</i>	<i>53</i>
<i>DROŠĪBAS NORĀDES</i>	<i>LV</i>	<i>55</i>
<i>NURODYMAI DĖL SAUGUMO</i>	<i>LT</i>	<i>56</i>
<i>УКАЗАНИЯ ЗА БЕЗОПАСНОСТ</i>	<i>BG</i>	<i>58</i>

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DE

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CE-DECLARATION OF CONFORMITY

EN

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EEC and 73/23/EEC and conforms with the standards or standardized documents DIN EN IEC 61010-1, -2, -010 und DIN EN IEC 61326-1.

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Warranty

You have purchased an original **IKA** laboratory machine which meets the highest engineering and quality standards.

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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	Page
CE - Declaration of conformity	2
Warranty	2
Safety instructions	16
Correct use	17
Unpacking	18
Commissioning	18
Connecting the pH-probe	18
Replacing extension cable and temperature sensor	18
Factory setting	19
Configuration	19
- Language	19
- Operating mode	19
- Maximum temperature	20
- Target temperature	20
- Temperature controller	20
- Time out	20
- pH-Measurement	21
- Reset	21
- Exit	21
- pH-Calibration	22
- PT 1000-Calibration	22
Error codes	23
Accessories	24
Technical Data	24

Safety instructions

For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be

accessed by everyone.

- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the media to be processed.
- **Caution!** The supply voltage (Limited Power Supply) for the **ETS-D6** must be provided using an IEC 61558 compliant transformer with double or reinforced insulation. Alternatively, laboratory equipment that meets the requirements of IEC 61010 or NEC Class 2 may be used.
- **Caution!** This equipment must only be used with media that have a flash point higher than the temperature safety limit set for the magnet agitator and/or heating system used. The safe temperature limit must always be set to at least 25 °C lower than the fire point of the media used.
- **Burn hazard!** Take care when touching the temperature sensor! The temperature sensor may heat up to 450 °C when inserted in the media!
- Do not touch the temperature sensor while measurements are being taken. This will prevent incorrect results.
- Ensure that the external temperature sensor (PT 1000, ETS-D ...) is inserted in the media to a depth of at least 20 mm, see fig. 2.
- Please ensure that the spiral cable does not come into contact with the heating plate.
- The stainless steel temperature sensor must not be used with aggressive media such as acids, caustic solutions or distilled water, due to the risk of corrosion. The glass sensor should be used in such cases.
- Only use glass encapsulated temperature sensors for electrolysis procedures.
- Always use the extension cable H 70 when the media being processed produces vapour. This ensures that the control unit does not come into contact with the vapour.
- Beware of hazards due to:
 - flammable materials
 - combustible media with a low boiling temperature

- **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- It may be possible for wear debris from rotating accessory parts to reach the material being processed.
- **Caution!** The user must check that all conditions have been met for the safe and correct operation of the equipment before operating the **ETS-D6** with a heater that does have a DIN 12878 compliant connector, and therefore does not provide cut-out protection. (The cut-out protection is a mechanism that ensures that the electrical heating circuit is switched off should the contact thermometer (**ETS-D6**) be broken or if the connector is removed.).

For protection of the equipment

- Protect the appliance and accessories from bumps and impacts.

Use of "H+P" magnetic stirring systems (no break fuse protection included)

When the supply current for the **ETS-D6** is interrupted, the electronics for the heatable magnetic agitator system will take over the temperature control. This means that the magnetic agitator system will heat up to the temperature set (e.g. 300 °C). Readjust the set-point temperature for the magnetic agitator system or switch the heating off.

Use of "Corning" magnetic stirrer PC-400, PC-420, PC-600, PC-620

When using the **ETS-D6** in conjunction with a "Corning" Magnetic

Agitator it is absolutely essential that the rotary "heat" dial on the magnetic agitator is set to the "O" (off) position. The **ETS-D6** will only control the set temperature when it is in this position.

Warning: The heating plate temperature limiter is not active in this operating mode.

If the rotary "heat" dial on the magnetic agitator is not in the "O" position then the magnetic agitator will heat up to the temperature set, regardless of the settings on the **ETS-D6**.

Warning: **ETS-D6** is inactive.

The safety advices given above represent the current status of the known risks. Nonetheless, the user must check equipment for correct and safe operation before operating any system comprising the **ETS-D6** and a magnetic stirring/heating system.

Correct use

- **Use**

- Precise measurement and control of temperature

The instrument can be connected to any magnetic agitator or heating plate that has a contact thermometer connection, provided this conforms with the requirements listed in the "Technical Data", see chapter "Technical Data".

These are amongst others **IKA** magnetic stirrers and Heidolph magnetic stirrers. You can connect Corning magnetic stirrers by the adapter AD-C1.

- pH-measurement

An optional pH-probe can be connected using the built in BNC connector.

When measuring temperature and pH at the same time, the glass sheathed temperature sensor should be used in order to avoid any interference with the pH measurement. To obtain useful results for pH measurements when using earthed heating plates (stainless steel, silumin) an electrically isolated container (glass beaker) must be used.

It is not possible to carry out pH measurements in metal containers.

The **ETS-D6** is not suitable for diagnostic measurements for medical purposes.

- **Range of use**

- Laboratories
- Pharmacies
- Schools

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.

Unpacking

- **Unpacking**

- Please unpack the device carefully
- In the case of any damage a detailed report must be set immediately (post, rail or forwarder)

- **Content of package**

- **ETS-D6**
- Operating instructions

Commissioning

☞ Plug the instrument into the contact thermometer jack on the heating magnetic stirrer or heating device , e.g.

IKA: RCT basic, IKA RET basic, RH basic KT/C,
RH digital KT/C, HCT basic

Heidolph: MR3001K, MR3002

VWR: VMS-D, VWR VMS-A

H+P Variomag

Please observe all safety instructions in this respect.

Connecting a "Corning" magnetic stirrer:

☞ Turn the rotary "heat" dial on the magnetic agitator to the "O" (off) position. Next, connect the **ETS-D6** to the magnetic agitator using the AD-C1 adapter provided.

Please observe all safety instructions in this respect.

➤ The device will be activated when the magnetic agitator is switched on.

If the supply voltage at the contact thermometer connector is not between 8 Vdc and 16 Vdc then it will be necessary to power the device using the H 52 power supply set (H 50 double connector adapter and H 51 plug-in power supply unit). The power supply and/or double connector adapter will provide the necessary supply voltage for the device. An automatic self-test will be performed when the device is activated. An introductory screen will be shown on the LCD display, see fig. 5.

Connecting the pH-probe

☞ Connect the BNC-connector for the pH-probe to the BNC-socket on the rear of the unit and lock in place using the bayonet fitting, see fig. 3.

Replacing extension cable and temperature sensor

☞ Pull the sensor and protective cap downwards, see fig. 4

➤ The plastic snap fastener is released

☞ Push the sensor or extension cable and protective cap over the fitting on the instrument

➤ The connection is re-established

➤ LCD display see fig. 28

Factory setting

The **ETS-D6** is supplied with the following factory settings:

Maximum temperature	450 °C
Mode	A
pH-measurement	disabled
Controller	PID
pH-Calibration	OK
PT 1000-Calibration	-
Extension cable	-
Time Out Er 5	3 min
Language	English

Configuration

The unit is configured using the Setup menu:

☞ Switch on the device

- LCD-Display:
1. Intro, see fig. 5
 2. Information about the operating mode/ Maximum temperature and the request "OK to enter Setup" (Duration: 2 seconds), see fig. 6

☞ Press the SET button during the two second period.

Use the ▲ and ▼ buttons to navigate through the menu and select the desired menu item.

The SET button is used to activate a selected menu item and confirm changes, see figs. 7, 8, 9.

The currently active settings are shown under the corresponding menu item.

Select the menu item "Exit" to leave the menu.

To accept and store changes select the "Yes" option when prompted with "Confirm changes" on the next screen.

➤ LCD display see fig. 10

• Language

The following languages can be selected under the menu item "Language":

- | | |
|-----------|------------|
| - English | - German |
| - French | - Spanish |
| - Chinese | - Japanese |
| - Korean | |

➤ LCD display see fig. 11

Accept and store the selected language:

☞ Press button SET

• Operating modes

The operating mode can be set under the menu item "Mode": The **ETS-D6** can be run in three different operating modes (A, B, C). The mode and maximum temperature currently set will be displayed for two seconds after power on, see fig. 6.

Mode A

The set-point temperature is always set to -10°C.

The maximum temperature can only be adjusted if Mode A is active.

☞ Use the ▲ and ▼ buttons to change the target temperature.

Mode B

The last target temperature selected will be used.

☞ Use the ▲ and ▼ buttons to change the target temperature

Mode C

The last values selected will be used.

☞ It is not possible to adjust the target temperature.

➤ LCD display see fig. 12

Accept and store the selected mode:

☞ Press button SET

Change the mode

☞ Hold down the ▲ button when turning on the unit. The system will switch to the next operating mode.

Sequence: A-B-C-A-B-C-A ...

Factory setting: Operating mode A

• **Maximum temperatur**

The maximum temperature (maximum allowed target temperature) can be set using the "Max. Temp" menu item.

☞ Use the ▲ and ▼ buttons to change the maximum temperature.
(0 - 450 °C)

➤ LCD display see fig. 13

Accept and store the selected maximum temperature:

☞ Press button SET

• **Target temperature**

The target temperature is set directly rather than using the Setup menu.

0,1 K steps

☞ Press buttons ▲ and ▼

1 K and 10 K steps respectively

☞ Press buttons ▲ and ▼ longer than 5 seconds

➤ LCD display see fig.14, 15:

- The instantaneous actual temperature with physical units
- The target temperature below the actual temperature
- The symbol  indicates heating is switched on
- The pH-value is shown at the bottom left if a pH-probe is connected

• **Temperature controller**

The type of control used can be set under menu item "Controller":

PID: Gives best results with minimum overshoot

On/Off: Operates with maximum heating speed

☞ Select Controller

➤ LCD display see fig. 16

Accept and store the selected controller:

☞ Press button SET

• **Time out**

"Time out" is the time taken to recognise that the sensor has not been placed in the medium when the heating is switched on.

Caution! This function will only be active if:

- Sensor temperature < 50 °C
- Difference target temperature/sensor temperature >5 K

The "Time out" period can be adjusted using the "Time-out Er5" menu item:

☞ Press buttons ▲ and ▼

➤ LCD display: "Time-out Er5" period see fig. 17

Accept and store the selected "Time-out" period:

☞ Press button SET

• pH-Measurement

The pH measuring function can be activated using the "pH-measurement" menu item.

➤ LCD display see fig. 18

Accept and store:

☞ Press button SET

• Reset (Restoring parameter to factory settings)

☞ Switch off the device

☞ Press and hold down the buttons ▲ and ▼, and switch on the equipment

or use the "Reset" menu item

➤ LCD display see fig. 19

Accept and store:

☞ Press button SET

• Exit

☞ Select the menu item "Exit" to leave the Setup menu.

➤ LCD display see fig. 9

Any changes made can be accepted or rejected.

pH-Calibration

The two-point calibration process can be started under menu item "ph calibration".

➤ LCD display see fig. 20

☞ Calibration can be cancelled at any time using the ▲ and ▼ buttons.

Calibration is used to adjust the unit for the measurement chain being used. The process determines the zero point and the pH gradient for the measurement chain.

The calibration should be carried out using pH 4 and pH 7 buffer solution as per DIN 19266.

Other buffer solutions, including those with identical rated values, may exhibit different temperature characteristics, which will lead to measurement errors.

- First, immerse the measurement chain and the temperature sensor in the pH 7 buffer solution. The rated value for the required buffer solution will be displayed, see fig. 20.
- Press the SET button to start the calibration process, see fig. 21 (duration up to 3 min.).
- Remove the measurement chain and the temperature sensor from the first buffer solution as soon as the rated value for the second buffer solution appears on display (see Fig. 22). Rinse both items thoroughly with distilled water and dry using a soft, lint free paper towel. Do not rub!
- Next, immerse the measurement chain and the temperature sensor in the pH 4 buffer solution.
- Press the SET button to restart the calibration process, see fig. 23.

The gradient for the measurement chain in % (mV/pH) will be displayed at the end of the calibration.

The calibration process was successful if the gradient is between 90%-105%, see fig. 24. The calibration values will be saved automatically.

- Confirm the setting using the SET button

If the gradient is outside this range then the calibration values will be rejected. The measurement chain should be swapped in this case.

PT 1000-Calibration

The **ETS-D6** is pre-calibrated without a sensor before shipping. Sensor tolerances in accordance with DIN IEC 751 Class A are not taken into account.

The complete measurement chain can be calibrated in order to eliminate the effects of the contact resistance at the plug connector and the sensor tolerances (inclusive of temperature drift).

The calibration process is started using the "PT1000 calibration" menu.

The following items are required:

- An additional temperature measuring device with an accuracy of + 0.05 °C
- Two stable, constant temperature sources (media at different, controlled temperatures)
- Immerse both sensors in the first medium. Wait until the display on the temperature measurement device has stabilised.
- Enter the measured temperature on the **ETS-D6** using the ▲ and ▼ buttons.
- Confirm the value entered using the SET button, see fig. 25.
- Immerse both sensors in the second medium. Wait until the display on the temperature measurement device has stabilised.
- Enter the measured temperature on the **ETS-D6** using the ▲ und ▼ buttons.
- Confirm the value entered using the SET button, see fig. 26.

The calibration is now successfully completed, see fig. 27.

Error codes

Any malfunctions during operation will be identified by an error message on the display, see e.g. fig. 29.

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- ☞ Disconnect the power supply for at least 10 sec.
- ☞ Carry out corrective measures
- ☞ Restart device

Error code	Cause	Effect	Correction
Er 1	Sensor not connected or defective (Sensor rupture)	Device off	- Connect temperature sensor - Replace temperature sensor
Er 3	Temperature inside device higher than 76 °C	Device off	- Switch off the device and allow to cool down
Er 5	Sensor not in media (no change in temperature of the media after heater active for 3 minutes)	Device off	- Insert the temperature sensor in the media to a depth of at least 20 mm. - Increase the "Time-out" period
Er 6	Temperature sensor short-circuit	Device off	- Switch off the device and replace temperature sensor
Er 7	The temperature of the medium is higher than the maximum target temperature that can be set	Device off	- Switch off the device and allow the media to cool down - Increase the maximum target temperature that can be set
Er 8	Calibration error	Device off	- Switch off the device - Re-calibration by authorised service personnel

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

Accessories

- H 38 Retaining bar
- H 16 Support rod
- H 44 Cross sleeve
- H 52 Power supply kit
- H 62.51 Standard sensor, ø 3 mm, length 270 mm, stainless steel
- H 70 Extension cable, length 1m
- H 66.51 Glass sensor, ø 8 mm, length 270 mm, glass
- AD-C1 Adapter (for "Corning" magnetic stirrer)

Technical Data

Temperature

Range of temperature measurement	°C	-50 - 450
Resolution	K	0,01
Measuring accuracy	K	± 0,05 + Fühlertoleranz PT 1000, DIN IEC 751 Klasse A, ± (0,15 + 0,002 x t), t in °C
Setting accuracy	K	0,1
Control deviation	K	± 0,2
<i>(500 ml water in 600 ml glass beaker, magnetic stirring bar 40, 600 min⁻¹, 50 °C)</i>		

pH-measurement

Range of measurement	pH	0 ... 14
Accuracy	pH	± 0,1
Resolution	pH	± 0,01
pH-connection		BNC jack
Supply voltage	Vdc	8 - 16 (Limited Power Supply)
Current consumption	mA	30 (bei 9V)
Perm. ambient temperature	°C	0 - 60
Perm. relative humidity	%	80
Perm. duration of operation	%	100
Protection type acc. to DIN EN 60529		IP54
Length spiral cable	mm	350
Male connector		6-pin DIN 45322
Socket		DIN 12878 Klasse 2
Dimensions B x T x H	mm	100x45x95 (without sensor)
Weight	kg	0,2
Operation at a terrestrial altitude	m	max. 2000

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