35 197 00 a

Dispergierantrieb T 10 basic Disperser T 10 basic

Betriebsanleitung	DE	4
Operating instructions	EN	11
Mode d'emploi	FR	18


Fig. 1


Fig. 2


Fig. 3

Pos.	Bezeichnung	Item	Désignation	Item	Désignation
1	Ein-/ Ausschalter	1	On-/ Off-switch	1	Interrupteur à deux positions
2	Drehzahlrad	2	Speed wheel	2	Roue de vitesse
3	Schaftverriegelung	3	Shaft locking	3	Verrouillage de la tige
4	Dispergierwerkzeug	4	Dispersion tool	4	Outils dispersantes
5	Klemme	5	Clamp	5	Collier
6	Drehgriff	6	Turning handle	6	Poignée
7	Kreuzmuffe	7	Cross sleeve	7	Manchon en croix
8	Stativ	8	Stand	8	Support

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DE

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Contents

	Page
CE - declaration of conformity	3
Warranty	11
Safety instructions	11
Correct use	13
Unpacking	13
Useful facts	13
Drive unit	14
Assembly of clamp and stand	14
Working with the unit	14
Error correction	15
Maintenance	15
Note to the tool set	15
Accessories	15
Permitted dispersing tools	16
Speed table	17
Technical data	17

Warranty

In accordance with IKA warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Safety instructions

In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

Only staff who have been trained accordingly, know the appliance and are authorised to carry out work in this field should use this appliance.

CAUTION! Never operate the equipment with the dispersing tools rotating freely. Ensure that parts of the body, hair or items of clothing cannot be trapped by the rotating parts. Before switching on the dispersing instrument make sure that its shaft is sufficiently immersed in the medium to prevent the medium from splashing out (see operating instructions of the dispersion tool). Please ensure, however, a minimum distance between the dispersing element and the base of the stirring vessel (fig. 1). Do not use the appliance without a dispersing element.

Vibrations can cause the turning handles on the clamp to loosen. Check that the turning handles are secure and tighten if necessary.

It is not permitted to use the dispersing instrument with a freely rotating dispersing element without a stirring vessel as any contact can cause serious injury. It is important to ensure that the dispersing element is always inside the stirring vessel when the appliance is switched on.

CAUTION! The equipment is also suitable for manual operation. Please follow the relevant safety instructions and guidelines, and occupational health and safety regulations for use in the laboratory. When working with the dispersing unit, the user must select and wear his personal protective equipment according to the mixing hazard category. Defective or inappropriate protective equipment can expose the user to the risk of splashing or getting caught in the dispersing tool.

Please note that the dispersing element and the journal bearings can become extremely hot during use.

Protect the dispersing element from any impact. Even small, undetectable damage creates unbalance and blocks the dispersing instrument. It also causes the dispersing element to overheat and inevitably overloads the drive.

Careful handling will guarantee safe work and a long service life of the machine.

In the event of unbalance or unusual noises, switch off the appliance immediately in order to prevent unwanted resonance both in the appliance and the whole assembly.

Any unbalance can damage the assembly, the appliance or destroy the stirring vessel, which, in turn, can seriously injure the user, as can the rotating dispersing element.

Unusual noises indicate that the appliance is no longer working properly.

Replace the dispersing element if you notice any unbalance or unusual noises. Stop using the appliance if there is no difference after you have changed the dispersing element and return it to the dealer or the manufacturer along with a description of the fault. Only dispersing and mixing elements approved by IKA may be used.

Even in the case of repair work, the unit must only be opened by a qualified engineer. The mains plug must be pulled before opening the unit. The voltage-carrying components inside the unit may continue to carry a voltage for some time after the mains plug is pulled.

With the hand enterprise and with the enterprise at the stand, it is to be made certain that safe conditions, a suitable anti-twist plate of the agitating container and a minimum distance of the dispersion tool are ensured to the agitating container soil.

Before switching on the dispersing instrument make sure that its shaft is immersed min. 20 mm in the medium to prevent the medium from splashing out.

The distance between the dispersion tool and the vessel bottom should not be less than 10 mm.

CAUTION! Covering or parts that are capable of being removed from the unit without accessory equipment have to be reattached to the unit for safe operation in order to prevent, for example, the ingress of fluids, foreign matter, etc..

Do not cover the ventilation slots on the drive in order to ensure adequate cooling of the drive.

The agitated vessels used for the dispersing have to be secured. Consider on a good stability of the entire structure.

Furthermore, the dispersing unit may only be used to stir those materials or material mixtures that the user knows will not react dangerously to the extra energy produced by the dispersing. This also applies to extra energy by means of solar radiation during the stirring procedure.

Do not process any pathogenic materials in open vessels as this may cause the user to inhale vapours. Always work under a suitable extractor hood, even when using closed vessels. Please contact IKA if you have any questions regarding work with pathogenic materials or other materials which require closed systems.

Make certain that the unit is set at the lowest speed before commissioning; otherwise, the unit will begin running at the speed last set. Dispersing at high speeds can cause the medium to splash out of the mixing vessel. There may be dangerous electrostatic activity between the medium and the output shaft.

The appliance starts up again following a cut in the power supply. Please note that the appliance must be switched off at the on/off switch in this case before the power supply is reconnected. The dispersing instrument may not be used in explosive rooms, for mixing dangerous substances or under water.

Please note that abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on. If there are any questions please contact IKA.

Always disconnect the plug before changing dispersing elements and attaching accessories.

The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug. The socket for the mains cord must be easily accessible.

Please observe the permitted speed for the dispersing element. Never set higher speeds.

Meter powder not too close to the flange, otherwise powder can be sucked into the drive.

Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.

Correct use

The dispersion equipment T 10 basic is suitable by use of most different dispersion tools for the production of emulsions or dispersions. It is conceived for the employment in laboratories.

The T 10 basic is conceived for the enterprise at a stand or for the enterprise at hand.

To the intended use the equipment must be fastened carefully over the arm with a cross sleeve to the stand. Make sure that the attachment mechanism (e.g. stand R 104) is set up nontilting on a even and sturdy underbody and during the dispersing procedure not to move begins.

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default. Do not touch rotating parts during operation!

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope includes:

One T 10 basic dispersion equipment according to the ordered type, one clamp R 200, a tool set and the operating instructions, all stored in a suitcase.

Useful facts

You have acquired a high- quality product. When used in combination with one of our recommended dispersing elements, the drive unit is a high-speed dispersing and emulsifying unit capable of handling free-flowing and liquid media in batches.

Dispersion is the dissolution and diffusion of a solid, liquid or gaseous phase in a continuum that is not consolute with that phase.

The rotor/stator principle:

Due to the high rotation speed of the rotor, the medium to be processed is automatically drawn axially into the dispersion head and then forced radially through the slots in the rotor/stator arrangement. The high accelerations acting on the material produce extremely strong shear and thrust forces. In addition, high turbulence occurs in the shear gap between rotor and stator, which provides optimum mixing of the suspension.

The dispersion effectiveness is heavily dependent on the product of the shear gradient and the time the particles spend in the shear zone. The optimum range for the circumferential velocity of the rotor/ stator arrangement is 6-20 m/s.

A processing time of a few minutes is usually sufficient to produce the desired fineness. Long processing times bring only insignificant improvements in the obtainable fineness; the energy expended serves merely to increase the temperature of the medium.

Drive unit

The open up a wide range of possibilities for dispersion technology in the conventional laboratory with a respective performance output of approximately 75 watts at 30.000 1/min.

The speed of the T 10 basic dispersion unit set by using the differently adjusting wheel.

Assembly of clamp and stand

Overview (fig. 1 and fig. 2)

Assembly

- Clamp (5) as in fig. 1 represented over the dispersion equipment push, until it engages in the position planned for it audibly.
- Tighten the turning handle (6) by hand up to the notice.
- Screw the clamp (5) (fig. 1) by means of cross sleeve onto the stand.

For safe working the drive units are fastened by means of a cross sleeve (7) to the plate stand (8). Pay attention when the assembling of the dispersion equipment at the stand to a perfect and firm seat of the cross sleeve.

In order to increase the stability of the mechanical structure, the drive unit must be installed as closely as possible to the stand seaweeds.

Disassembly

- For the disassembly of the dispersion equipment of the clamp (5), first the turning handle (6) is to be loosened.
- Hold the dispersion equipment with a hand, with which other hand now the dispersion equipment knows through pressures on the locking lever (fig. 2) from the clamp is solved.

Working with the unit

Overview (fig. 1)

Check that the voltage listed on the designation plate matches the mains voltage. Also pay attention to the ambient conditions listed in the technical data. The power socket used must be earthed (protective earth conductor contact).

Assembly of the dispersion tool

- Open the shaft lock (3) by turning to the left as far as the stop.
- Insert the dispersing element into the collet as far as the stop (you will feel it click into place after some resistance).
- Lock the shaft lock (3) by turning it to the right until it noticeably clicks into place. The marks on the shaft lock (3) and on the casing should be over each other when closed.

The distance between the dispersion tool and the vessel bottom should not be less than 10 mm.

Please observe also the operating instructions of the dispersion tool to ensure a safe usage.

The entire unit may also be arranged somewhat off centre in order to avoid any undesired air admission caused by the strong rotational turbulence.

Set the adjusting wheel (2) to the minimum speed before switching on the appliance. You can read off the set speed by comparing the numbers on the adjusting wheel (2) and the speed table (see chapter: Speed Table).

By using the adjustment wheel (2) the speed may be infinitely varied to meet the requirements of the medium to be processed. The no-load speed of the T 10 basic driving unit may be set between 8.000 and 30.000 1/min.

The unit is ready for service when the mains plug has been plugged in.

The driving unit is switched on by means of the rocker switch (1).

Error correction

The equipment is appropriate for short-time duty (10 min ON / 5 min OFF), i.e. indicated fall times must be considered. During neglect the equipment can switch off depending upon load by temperature rise during the enterprise, since the drive is equipped with an overload safety device, which switches the equipment off with overload or temperature rise. If the equipment switched off at the On-/ Off- switch. After an appropriate cooling time the equipment can be restarted.


Note: The equipment does not restart it after responding the overload safety device automatically, must for cooling at the rocker switch (1) be switched off.

Maintenance

Maintenance and Cleaning

Drive unit:

The drive unit is maintenance-free but not immune from wear. The carbon brushes of the motor wear down over the time.

For cleaning disconnect the main plug

To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Dispersion tools:

Consider please the manuals associated to the dispersion tools.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see **www.ika.net**, spare parts diagram and spare parts list

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Note to the tool set

The provided assembly tools (flat key and shaft key, fig. 3) can be used to disassemble and assemble the dispersing elements for cleaning and maintenance. Take the information about the exactly use of the operating instructions of the dispersion tools.

Accessories

- | | |
|----------|--------------|
| • R 1825 | Plate stand |
| • R 104 | Stand |
| • R 182 | Cross sleeve |
| • H 44 | Cross sleeve |

Permitted dispersion tools

Designation		Material of shank
S10N-5G		Stainless steel
S10N-8G		Stainless steel
S10N-10G		Stainless steel
S10D-7G-KS-65		Synthetic material
S10D-7G-KS-110		Synthetic material

Only use the dispersion tools listed in the table and note the according operating instructions of the dispersion tool.

Speed table

Note: The speed table was calculated on the basis of a dispersion volume of 5 ml water using the dispersing tool S10N-5G.

With other combinations (dispersion tools; medium quantity and medium viscosity) can result other number of revolutions values.

Adjusting wheel scale	1	2	3	4	5	6
No-load speed	8.000	9.500	11.500	14.500	20.500	30.000
Speed with 5 ml water	7.900	9.400	11.400	14.450	20.450	29.900

The speeds depend on the dispersing element used as well as the viscosity and quantity of the medium. It may, therefore, change over the course of a dispersing process if the viscosity of the liquid changes.

Technical data

Speed range	rpm	8.000 - 30.000 (at nominal design voltage and 50 Hz)
Speed display		Scale (see speed table)
Speed variation on load change	%	<6
Perm. ambient temperature	°C	5 - 40
Perm. relative humidity	%	80
Perm. on time (drive unit)	min	max 10 ON min. 5 OFF
Overload protection / blocking protected		locking bimetallic switch; temperature and power sensitively
Protection to DIN EN 60529		IP30
Power consumption	W	125
Power output	W	75
Voltage	VAC	220 - 240 (nominal design voltage 230V)
	or	VAC 115
	or	VAC 100
Frequency	Hz	50/60
Noise (without dispersing element)	dB(A)	65
Dimension (W x D x H) drive	mm	46x57x201
Dimension clamp arm	mm	Ø8 / L100
Weight	kg	0,4
Contamination level		2
Protection class		II
Overvoltage categorie		II
Operation at a terrestrial altitude	m	max. 2000 above sea level

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GmbH & Co.KG**
IKA® Works Inc.