IKA® REO basic C

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CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 61 010; EN 50 082; EN 55 014; EN 55 114 and EN 60 555.

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Staufen, Mai 2002

Reiner Dietsche
Vice President

Wolfgang Buchmann
Corporate Director - Quality Assurance

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Garantee

You have purchased an laboratory machine which meets the highest engineering and quality standards.

In accordance with guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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Contents

| | Page |
|---|------|
| Guarantee | 2 |
| Safety instructions | 6 |
| Proper use | 6 |
| Unpacking | 6 |
| Commissioning | 7 |
| Switching on | 7 |
| Function Stirring | 7 |
| Maintenance and Cleaning | 7 |
| Associated standards and regulations | 8 |
| Accessories | 8 |
| Technical data | 8 |
| Information for care and maintenance of the heating plate vers. C | 12 |
| List of spare parts REO basic C | 14 |
| Spare parts diagram REO basic C | 15 |

Safety instructions

The instrument may not be operated in rooms with explosion hazard.

Attention - Magnetism! Effects of the magnetic field have to be taken into account (e.g. data carriers, cardiac pacemakers ...).

When using PTFE-coated magnetic bars, the following have to be noted: *Chemical reactions of PTFE occur in contact with molten or dissolved alkaline and alkaline-earth metals, as well as with fine-particled powders of metals of the 2. and 3. group of the periodical system at temperatures above 300-400°C.*

Only elementary fluorine, chlorine trifluoride and alkaline metals do attack PTFE, halogen hydrocarbons have a reversibly swelling effect.

Source: Römpps Chemie-Lexikon and „Ullmann“ Bd.19

Even in the case of repair work, the unit must only be opened by a qualified engineer. The mains plug must be pulled before opening the unit. The voltage-carrying components inside the unit may continue to carry a voltage for some time after the mains plug is pulled.

Proper use

The IKAMAG REO basic C is a laboratory device. It is used in laboratories in the chemical industry, schools and pharmacies to name a few. The devices are suitable for mixing of substances which are placed on mounting plates. The mixing intensity depends on the motor speed and the size of the magnetic stirring bar.

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers: an REO C, connecting cable and operating instructions.

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.

Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

Switching on

After the instrument has been switched on, it is in „standby“-operation. The function „stirring“ is started with the help of switch on the right bottom of the housing. (see illustration 1)


illustration 1

Function Stirring

The funktion stirring is started and stopped by its corresponding switch. If the function “Stirring” is started, the LED “Mot” beside the the imprint “Mot” shines green.

The speed of the drive can be adjusted on the “Mot” dial with the corresponding scale from 50 to 1500 rpms.

The motor speed of the instrument is kept constant by a control circuit.

Maintenance and Cleaning

The REO basic C is maintenace-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

| | |
|------------------------|--|
| Dyes | isopropyl alcohol |
| Construction materials | water containing tenside / isopropyl alcohol |
| Cosmetics | water containing tenside / isopropyl alcohol |
| Foodstuffs | water containing tenside |
| Fuels | water containing tenside |

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

Protective conductor:

On the under side of the mounting plate are areas (bare metal) and can be used as a test point for the protective conductor. It may be necessary to remove an oxide coating for the subsequent tests.

Associated standards and regulations

Construction in accordance with the following safety standards

EN 61 010-1

EN 61 010-2-10

UL 3101-1

CAN/CSA C22.2 (1010-1)

Construction in accordance with the following EMC standards

EN 61 326-1+A1

Associated EU guidelines

EMC-guidelines: 89/336/EWG

Machine guidelines: 73/023/EWG

Accessories

RS 1 set of stirring bars **H16V** support rod

RSE PTFE stirring bar remover

H 44 cross sleeve

H 36 holding rod

Recommended IKA-stirrers

TRIKA-stirrers PTFE, triangular
lengths 25, 40 and 55 mm

IKAFLON-stirrers PTFE

dia 6mm lengths 10 and 15 mm

dia 8mm lengths 20, 25, 30, 40, 50 and 60 mm

Technical data

| | |
|---|---|
| design voltage: | VAC 220 - 240 ±10% |
| | oder VAC 115 ±10% |
| | oder VAC 100 ±10% |
| design frequency: | Hz 50/60 |
| Drive: | EC-motor with ball bearings, without brushes |
| Max. power consumption / output: | W 12 / 5 |
| speed range: | rpm 50...1500 infinitely variable setting |
| speed variety (without load, 1500rpm): | % -5....+15 |
| stirring quantity (H ₂ O): | l 20 |
| diameter of mounting plate: | mm 135 |
| permitted ambient temperature: | °C +5 to +40 |
| permitted relative humidity: | % 80 |
| permitted duration of operation: | % 100 |
| protection class acc. DIN 40 050: | IP 42 |
| fuses: | F1/F2 10AT/250V Id-Nr 14 080 00 |
| contamination level: | 2 |
| overvoltage category: | II |
| protection class: | 1(protective earth) |
| operation at a terrestrial altitude | m max 2000 above sea level |
| dimensions: | mm 160 x 280 x 85 |
| weight: | kg 2,4 |

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Information for Care and Maintenance of the Mounting Plate with Technical Enamel Coating


The technical enamel coating makes the mounting plate easier to care for and more resistant to acids and bases. Because of it, however, the mounting plate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.

Make certain that the bottom of the placing vessel is even, clean and dry. The bottom of the placing vessel must not have any sharp grooves, sides or edges. Remove residues of bases and acids immediately. We recommend most strongly that you clean the mounting plate regularly.

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
List of spare parts REO basic C

| Nr. | Designation |
|------|---------------------------|
| 10 | housing-upper part |
| 11 | potentiometer shaft |
| 12 | clamping bushing |
| 13 | spring nut |
| 20 | washer |
| 30 | EC - motor |
| 31 | housing-bottom part |
| 32 | plastite screw |
| 33 | PCB control |
| 34 | bushing |
| 35 | magnet |
| 36 | countersunk screw |
| 37 | plastic cover left |
| 38 | taptite countersunk screw |
| 39 | cap |
| 40 | screwed plug |
| 41 | operation knob |
| 44 | cover |
| 45 | taptite chese-head screw |
| 46 | plastic cover right |
| 47 | fixing plate |
| 50 | lock washer |
| 101 | mounting plate |
| 108 | spring nut |
| 3001 | poti 1k |
| 3002 | electric mains |
| 3003 | rocker switch motor |
| 3007 | trafo |
| 3008 | fuse 6,3AT |
| 3010 | connector 3pol |
| 3011 | connector 4pol |
| 5001 | grooved ball bearing |
| 5002 | grooved ball bearing |

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IKA® WORKS, INC.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

The Global Market of Laboratory Equipment

<http://www.labworld-online.com>

Europa - Afrika

IKA® -WERKE GMBH & CO.KG

Janke & Kunkel-Str. 10
D 79 219 STAUFEN
GERMANY
TEL. 07633/831-0
FAX 07633/831-98
E-mail: sales@ika.de
Internet: <http://www.ika.net>

America

IKA® WORKS, INC.

2635 NORTH CHASE PKWY. SE
WILMINGTON, NC 28405-7419
TEL. 800/733-3037
TEL. 910/452-7059
FAX 910/452-7693
E-mail: usa@ika.net

Asien - Australien

IKA® Works (Asia) Sdn Bhd

(Company No. 340448-K)
Lot 2, Jalan Indah 1/2
Taman Industri Rawang Indah
48000 Rawang
Selangor, Malaysia
TEL :(603) 6093 3322
FAX :(603) 6093 3940
E-mail: ika@tm.net.my