

## Maxi MR 1 digital IKAMAG®



- Powerful magnetic stirrer without heating.
- Flat, sturdy stainless steel casing
  - Non-locking motor
  - Infinitely variable speed
  - Digital LED speed display
  - Timer (0 - 56 min) or continuous operation
  - For stirring quantities up to 150 liters (H<sub>2</sub>O)

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H <sub>2</sub> O) [l]	150
Motor rating input [W]	80	Motor rating output [W]	35
Speed display [-]	LCD	Speed range [1/min]	0 - 600
Stirring bar length max. [mm]	155	Set-up plate material [-]	stainless steel 1.4301
Set-up plate dimensions [mm]	500 x 500	Dimensions (W x H x D) [mm]	505 x 110 x 585
Weight [kg]	16	Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 21
Voltage [V]	230 / 100 - 120 / 100	Frequency [Hz]	50/60
Power input [W]	80		