

## RO 5 power IKAMAG®



Multi-position magnetic stirrer with 5 stirring positions, without heating. The stainless steel surface covers the unit allowing easy cleaning and providing protection against the penetration of liquids.

- Optimum use of laboratory space
- Including removable PUR cover

<b>Technical Data</b>			
Number of stirring positions [-]	5	Stirring quantity max. per stirring position (H <sub>2</sub> O) [l]	0.4
Stirring position distance [mm]	90	Speed deviation of individual stirring positions [%]	5
Motor rating input [W]	7.2	Motor rating output [W]	1.8
Speed range [1/min]	0 - 1100	Stirring bar length max. [mm]	30
Speed control [-]	scale 1-10	Set-up plate material [-]	stainless steel 1.4301
Set-up plate dimensions [mm]	120 x 450	Dimensions (W x H x D) [mm]	122 x 65 x 552
Weight [kg]	2.3	Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 42
Voltage [V]	230 / 115	Frequency [Hz]	50/60
Power input [W]	10		