

KS 130 basic



Small, quiet shaker ensures long life with ideal swivel motion, for a maximum shaking weight of 2 kg.

- Electronic adjustment of speed and timer
- LED display for speed and time adjustment
- Wide range of attachment combinations makes it possible to use almost all shapes and sizes of vessels
- Attachments are not included, please order separately

| Technical Data | | | |
|--|--------------------------------|-----------------------------------|-----------------------|
| Type of movement [-] | orbital | Orbital diameter [mm] | 4 |
| Permissible shaking weight (incl. attachment) [kg] | 2.0 | Motor rating input [W] | 45 |
| Motor rating output [W] | 10 | Permissible ON time [%] | 100 |
| Speed min (adjustable) [1/min] | 80 | Speed range [1/min] | 0 - 800 |
| Speed display [-] | LED line | Timer [-] | yes |
| Timer display [-] | LED line | Time setting range [min] | 5 - 50 |
| Operating mode [-] | timer and continuous operation | Shaking table lock [-] | no |
| Dimensions (W x H x D) [mm] | 270 x 98 x 316 | Weight [kg] | 8.8 |
| Permissible ambient temperature [°C] | 5 - 50 | Permissible relative moisture [%] | 80 |
| Protection class according to DIN EN 60529 [-] | IP 21 | RS 232 interface [-] | no |
| Analog output [-] | no | Voltage [V] | 220 - 240 / 100 - 120 |
| Frequency [Hz] | 50/60 | Power input [W] | 45 |